





Organic hydride



Organic hydride



Organic hydride

Conditions

System	: GC-FID
Column	: InertCap AQUATIC-2 0.32 mm I.D. x 30 m df = 1.8 μ m
Col. Cat. No.	: 1010-19247
Col. Temp.	: 40 °C (3 min hold) – 5 °C/min – 90 °C – 10 °C/min – 220 °C
Carrier Gas	: He 1.7 mL/min
Injection	: Split flow 34 mL/min 230 °C
Detection	: FID Auto Range 250 °C
Sample Size	: 0.5 μ L