


Phenothiazins



Column : Inertsil ODS-3
 5µm
 150 × 4.6mm I.D.

Column Cat. No. : 5020-13001

Eluent : CH₃CN/0.1M KH₂PO₄(pH 2.9 H₃PO₄)=1/2

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 254 nm

Injection Vol. :

Sample :

Solutes : 1. Prochlorperazine
 2. Chlorpromazine
 3. Fluphenazine
 4. Triflupromazine
 5. Thioridazine