


Analysis of Fluocinolone acetonide(column select)

Data No. LA171-0000



$R_s=2.5$

Column :Inertsil ODS-3V
5 μ m
250 \times 4.6mmI.D.

Solutes : 1. iso-Propyl paraben(0.05mg/mL)
2. n-Propyl paraben(0.05mg/mL)

Column Cat. No. :

Eluent :CH₃CN/H₂O=30/70

Flow rate :2.0 mL/min.

Col. Temp. :40°C

Detection :UV 254 nm

Injection Vol. :20 μ L

Sample : 1. iso-Propyl paraben(0.05mg/mL)
2. n-Propyl paraben(0.05mg/mL)