


Polyaromatic hydrocarbons



Column : Inertsil ODS-2
 5µm
 250 × 4.6mm I.D.

Column Cat. No. : 5020-01102

Eluent : CH₃CN/H₂O=75/25

Flow rate : 1.0 mL/min.

Col. Temp. : 40°C

Detection : UV 220 nm

Injection Vol. :

Sample :

Solutes : 1. Naphthalene
 2. Acenaphthylene
 3. Biphenyl (Diphenyl)
 4. Fluorene
 5. Acenaphthene
 6. Phenanthrene
 7. Anthracene
 8. Benzo(a)anthracene
 9. Fluoranthene
 10. Pyrene
 11. Triphenylene
 12. Crycene
 13. Benzo(e)pyrene
 14. Benzo(a)pyrene
 15. Dibenzo(ah)anthracene
 16. Benzo(ghi)perylene