


Analysis of aliphatic acid (prelabel derivatization by ADAM reagent)



Conditions

System	: GL-7400 HPLC system
Column	: Inertsil ODS-P (5 μ m, 250 x 3.0 mm I.D)
Column Cat. No.	:
Eluent	: CH ₃ CN
Flow Rate	: 0.9 mL/min
Col. Temp.	: 20 °C
Detection	: FL Ex 365 nm Em 412 nm (GL-7453A FL detector)
Injection Vol.	: 5 μ L
Analyte	: aliphatic acid derivitized with ADAM reagent