


Analysis of Carbendazim (Inertsil Ph-3)

**Conditions**

System : GL-7400 HPLC system
Column : Inertsil Ph-3 (5 μ m, 150 x 4.6 mm I.D.)
Column Cat. No. : 5020-01920
Eluent : A) CH₃CN
 : B) H₂O
 : A/B = 25/75, v/v
Flow Rate : 0.7 mL/min
Col. Temp. : 40 °C
Detection : FL Ex 285 nm Em 315 nm (GL-7453A FL Detector)

Analyte

1. Carbendazim (1 mg/L)